





TECHNICAL DATA SHEET

FORMALDEHYDE CHEMICAL INDICATOR

TECHNICAL DATA SHEET	PRODUCT CODE	:	4AFOI	
	RELATED STANDARDS	:	4A Medical Formaldehyde Chemical Indicator is classified as Type 4 according to ISO 11140-1 standard.	
	INTENDED USE	:	It is indicator for package which is used in formaldehyde sterilization, it is produced according to ISO 11140 Standards. Indicators which are included in Type 4 according to ISO 11140 standards gives information to user whether material within package is sterile or not.	
	OPERATION CONDITIONS	:	It is used to evaluate the process of formaldehyde sterilization. 4A Medical Formaldehyde indicators are water based and they do not include toxic substances.	
			<u>Before Sterilization</u>	<u>After Sterilization</u>
				
				
	PRODUCT SPECIFICATIONS	:	Each box contains 500 pieces indicator printed carton sheets. The indicator is formaldehyded concentration sensitive.	
	STORAGE AND DISPOSAL	:	The expired products are not used. It's kept in the places with room temperature without humidity and sun lights. It shouldn't be in contact with liquid. Damaged and defective products are not used. Shelf life is 5 years after production date. After sterilization indicator keep sterilization colour for 6 month.	
	PACKAGING	:	20 packs in cardboard.	